


NNIN Nanotechnology Education

How Does a Light Compound Microscope Work?

A compound light microscope contains two sets of lens which increases magnification. Normally light bounces off an object in a straight line. In a microscope the lens causes the light waves to bend in toward each other forming a "cone" of light which focuses on the next lens. When the light reaches the eye the object has been magnified to appear hundreds of times larger than its original size.





Images from: http://www.cas.miamioh.edu/mbi-ws/microscopes/compoundscope.html

www.nnin.org

National Nanotechnology Infrastructure Network

Copyright University of xxxx 200x

Permission granted for printing and copying for local classroom use without modification

Development and distribution partially funded by the National Science Foundation

NNIN Document: NNIN-1336

Rev: 08/2013